

ICE LINED REFRIGERATOR I525H120



**ENERGY
CONSUMPTION**
0.68 kWh/24 hrs.

SPECIFICATIONS

Gross Volume:	266 L
Vaccine Capacity:	216 L
Dimension (mm) WxDxH:	1410 x 805 x 970
Temperature Range:	2°C to 8°C
Operating Condition:	10°C to 43°C
Refrigerant:	R-290

FEATURES



Copper
Evaporator Coil



SS304
Internal Cabinet
for Rust-free Life



Controller with
72hr Display backup,
High-temp Alarm



High ICE Pack
Freezing Capacity



Long Hold Over
Time with 80mm
PU foam walls



Power cord
with Anti-rodent
Conduit pipes

