

PHARMA REFRIGERATOR PR600



ENERGY
CONSUMPTION
1.23 kWh/24 hrs.

SPECIFICATIONS

Gross Volume:	507 L
Vaccine Capacity:	453 L
Dimension (mm) WxDxH:	660 x 830 x 1817
Temperature Range:	2°C to 8°C
Operating Condition:	5°C to 43°C
No. of Shelves:	4
Refrigerant:	R-290

FEATURES



SS304
Internal Cabinet
for Rust-free Life



Digital
Temperature
Display Indicator



Power cord
with Anti-rodent
Conduit pipes



Maintenance
Free Condenser

