

VACCINE REFRIGERATOR VFW140H-HC



**ENERGY
CONSUMPTION**
1.94 kWh/24 hrs.

SPECIFICATIONS

Gross Volume:	142 L
Vaccine Capacity:	68 L
Dimension (mm) WxDxH:	710 x 720 x 955
Temperature Range:	-15°C to -25°C
Operating Condition:	5°C to 43°C
No. of Baskets:	4
Refrigerant:	R-290

FEATURES



Copper
Evaporator Coil



SS304
Internal Cabinet
for Rust-free Life



Controller with
72hr Display backup,
High-temp Alarm



High ICE Pack
Freezing Capacity



Long Hold Over
Time with 80mm
PU foam walls



Power cord
with Anti-rodent
Conduit pipes

